


B - LAMINA[®]

Composite double sided implants



Specifications

- Density of polypropylene (50g/m²) (parietal side)
- Silicone (visceral side)

Indications

B - Lamina[®] has been designed to be used through laparoscopic or conventional open surgery approach for hernia repairs, eventrations and as a reinforcement of the abdominal wall in intraperitoneal position.


Description

B - Lamina[®] is a non-woven 2D implant composed of a microporous structure made from 100% polypropylene fibers thermally bonded (parietal side) and silicone (visceral side).

Benefits

- 100% Biocompatible
- The silicone coating prevents the formation of adhesions in contact with the visceral peritoneum
- The microporous structure made from polypropylene fibers placed on the parietal side, provides excellent drainage and avoids the formation of "seroma".

Code	Sizes	Box
EMS12	Ø 12 cm	1
EMS1015	10 x 15 cm	1
EMS1517	15 x 17 cm	1
EMS1522	15 x 22 x 10 cm	1
EMS2030	20 x 30 cm	1
EMS3030	30 x 30 cm	1





Italia

Stabilimento di Produzione

Zona Industriale
67062 Magliano dei Marsi (AQ)
Tel. +39 0863 515215/6 - Fax +39 0863 570070

Sede Legale

Via G. Gregoraci, n.12
00173 Roma
Tel. +39 06 72677319 - Fax +39 06 72679909

France

86-114 Avenue Louis Roche - Batiment D Hall-3
92230 Gennevilliers
Tel. +33 147985294 - Fax +33 147985892

España

Polgono Ind. Mutilva Baja - Calle V, n.16
31192 - Mutilva Baja
Tel. +34 948122274 - Fax +34 948120842

UK

Unit 12a Barncliffe Mills
HD8 8LU Huddersfield
Tel. +44 1484602222 - Fax +44 1484602280

USA

6600WE, Rogers Circle
Suite 1&2
Boca Raton, FL 33487

Republic of Belarus

V. Gostinets Street , 143
222310 - Molodechno - MINSK
Tel. +375 1767 45555 - Fax +375 296 374737

Brasil

Rua Assunção n. 411 Salas 101 e 201
Botafogo Rio de Janeiro (RJ)
CEP 22251 - 030
Tel. +5521 22661943 - Fax +5521 22663539

Algerie

Zone D'Activite dar el Beida
Lot n° 66, Bt 12
16100 Alger

Tunisie

Rue Ahmed Echouk El Menzah 16, 9a
1082 Tunis
Tel. +216 71336017 - Fax +216 71254912

www.ergonsutramed.com
info@ergonsutramed.com



CE 0373 - UNI EN ISO 9001:2008 - UNI EN ISO 13485:2003

Il presente documento è destinato esclusivamente ad uso degli operatori sanitari • Fare sempre riferimento al foglietto illustrativo e leggere attentamente le istruzioni per l'uso • I prodotti potrebbero non essere disponibili in tutti i mercati • Contatta un promotore Ergon Sutramed per chiarimenti o domande circa la disponibilità dei prodotti Ergon Sutramed.
This document is intended solely for the use of healthcare professionals • Please refer always to the instruction for use • Products may not be available in all markets • Please contact your Ergon Sutramed representative if you have questions about the availability of products in your area.